QL[®] Locking Syringe

Locking Syringe Device

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
Atrion Medical Products, Inc.
1426 Curt Francis Road
Arab, Alabama 35016 USA

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QL® Locking Syringe

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in serious complications.

Description:

The QL® Locking Syringe is a one-piece, plastic, disposable syringe device with a lock lever design that controls the piston and a connecting tube with a male rotating adapter.

Indications:

The QL® Locking Syringe is intended for single use to inject fluid into and aspirate fluid from a balloon interventional device having a recommended inflation volume.

Contraindications:

None

Warnings:

- Use only liquid inflation media. Do not inflate with air or gas.
- Always follow the manufacturer's directions accompanying the balloon catheter for instructions for use, precautions, and warnings for that device.

Precaution:

- Before use, inspect the device to verify that no damage has occurred during shipping and handling.
- Before use, ensure the fluid path, including the connector tubing is completely free of air.

Instructions for Use:

Preparation:

Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked.

Unlock the piston by pushing the lock lever left. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

1. Prepare a solution of contrast media and saline in a small sterile bowl. Check balloon interventional device and contrast media instructions for specific contrast mixture recommendations.
2. Orient the tubing downward into the contrast media.
3. Aspirate enough solution to fill the syringe.
4. Hold the device upright to purge the air from the syringe and connecting tube. Tap the syringe lightly, if necessary, to remove all the air bubbles and fill the connecting tube completely.
5. Inspect the syringe and tubing to ensure that the device has been completely purged of air.
6. Adjust the syringe volume recommended for the balloon interventional device. If more contrast solution is needed, submerge the syringe hose into the basin of solution and aspirate.

Attaching the inflation device to the balloon catheter:

1. Prepare and test the balloon catheter according to the manufacturer's directions for use.
2. If a separate syringe was used to prepare the balloon catheter, remove it. When a high pressure stopcock is installed on the end of the inflation device connecting tube, it should be opened and purged with contrast media from the inflation device to eliminate air. Create a fluid-fluid connection between the balloon and the stopcock or connecting tube (male rotating adapter) of the locking syringe by placing a drop of contrast solution from the syringe into each hub.
3. Hand-tighten the hubs securely.

Operating locking syringe:

1. Fill device with the recommended amount of solution, then slide the lock lever to locked position by pushing the lever to the right for fine volume adjustments and to maintain the piston at the desired volume setting.
2. To inject fluid into the balloon, release the lock lever and manually push the plunger forward until you reach the balloon manufacturer's recommended volume. Note: For slow injection, with the lock lever in locked position, turn the palm grip on the piston clockwise until the desired volume of fluid has been delivered.
3. For rapid aspiration, with the lock lever in its unlocked position, pull back the palm grip manually. For slow aspiration, with the lock lever in the locked position, turn the palm grip on the piston counterclockwise until the desired aspiration volume has been reached. The lock lever may be slid back to locked position to maintain the piston in the aspiration position.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes. Sterile unless package opened or damaged.

US Patent No. 5,713,242, 6,796,959B2, and 6,938,319B2

⚡ **Atrion** and QL® are Registered Trademarks of Atrion Medical Products, Inc.

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