QL[®] Fluid Dispensing Syringe

Fluid Dispensing Syringe

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Atrion Medical Products, Inc.
1426 Curt Francis Road
Arab, Alabama 35016 USA

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Manufacturer
製造者

QL® Fluid Dispensing Syringe

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in serious complications.

Description:

The QL® Fluid Dispensing Syringe is a one-piece, plastic, disposable syringe with a lock lever design that controls the piston, a manometer, and a connecting tube with a male rotating adapter. An optional 3-way stopcock may also be enclosed for use during preparation of the device. The manometer measures pressures ranging from vacuum to gauge capacity. The accuracy of the manometer has been determined to be within 5 psi over the range.

Indications:

The syringe is recommended for use by healthcare professionals to dispense fluids to the body and monitor the pressure of those fluids.

Contraindications:

None

Warnings:

- Use only liquid media. Do not inflate with air.
- Always follow the manufacturer's directions accompanying the other instruments needed for the relevant procedure in which fluids will be dispensed into the body and the pressure of those fluids monitored.

Precaution:

- Before use, inspect the device to verify that no damage has occurred during shipping and handling.
- Before use, ensure the connector tubing is completely free of air.

Instructions for Use:

Preparation:

Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked.

Unlock the piston by pushing the lock lever left. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

1. Prepare a solution in a small sterile bowl.
2. Attach an open stopcock to the end of the tubing and immerse downward into the solution.
3. Push the release lever left and aspirate enough solution to fill the syringe one-half to two-thirds full, then close the stopcock.
4. While holding the device with the distal end facing down, pull back plunger knob to create a vacuum and purge air from the gauge, syringe and connecting tube. Tap the syringe lightly, if necessary, to loosen all the air bubbles, then position device upright, open stopcock and push plunger forward to expel air and fill the connecting tube completely with fluid.
5. Inspect the syringe and tubing and stopcock to ensure that the device has been completely purged of air bubbles.
6. Adjust the syringe volume to the desired amount. If more solution is needed, submerge the tubing tip into the basin of solution and aspirate. Close stopcock.

Operating Syringe:

1. Release the lock lever and allow the piston to move forward into neutral position (0 psi).
2. To increase pressure, engage the lock lever, turn the palm grip on the piston clockwise slowly until the desired pressure is reached. The lock lever maintains the increasing pressure.
3. To gradually decrease pressure, turn the palm grip on the piston counterclockwise slowly until the desired pressure is reached.
4. To rapidly decrease pressure, push the lock lever left, releasing the piston, and pull back. Slide the lock lever back to lock, if desired.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

US Patent No. 5,713,242, 6,796,959B2, and 6,938,319B2

Atrion and QL® are Registered Trademarks of Atrion Medical Products, Inc.

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