QL[®] Fluid Dispensing Syringe

Fluid Dispensing Syringe

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QL® Fluid Dispensing Syringe

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in serious complications.

Description:

The QL® Fluid Dispensing Syringe is a one-piece, plastic, disposable syringe with a thumb ring design that controls the piston, a manometer, and a connecting tube with a male rotating adapter. An optional 3-way stopcock may also be enclosed for use during preparation of the device. The manometer measures pressures ranging from vacuum to gauge capacity. The accuracy of the manometer has been determined to be within 10 psi over the range.

Indications:

The syringe is recommended for use by healthcare professionals to dispense fluids to the body and monitor the pressure of those fluids.

Contraindications:

None

Warnings:

- Use only liquid media. Do not inflate with air.
- Always follow the manufacturer's directions accompanying the other instruments needed for the relevant procedure in which fluids will be dispensed into the body and the pressure of those fluids monitored.

Precaution:

- Before use, inspect the device to verify that no damage has occurred during shipping and handling.
- Before use, ensure the connector tubing is completely free of air.

Instructions for Use:

Preparation:

You can freely pull the piston back for aspiration, or push it forward for injection.

1. Prepare the solution in a sterile vessel according to the solution manufacturer's instructions.
2. Attach an open stopcock to the end of the tubing and bring it into contact with the desired solution.
3. Pull back on the plunger thumb ring to aspirate the desired volume of solution into the syringe; then close the stopcock.
4. While holding the device with the hose end facing down, pull back on the plunger thumb ring to create a vacuum and purge air from the gauge, syringe and connecting tube. Tap the syringe lightly, if necessary, to loosen all the air bubbles, then position device upright, open stopcock and push plunger thumb ring forward to expel air and fill the connecting tube completely with fluid.
5. Inspect the syringe and tubing and stopcock to ensure that the device has been completely purged of air bubbles.
6. Adjust the syringe volume to the desired amount. If more solution is needed, repeat the previous filling and purging procedures. Close stopcock.

Operating Syringe:

1. To deliver fluid, increase pressure or hold a set pressure, push plunger thumb ring forward until the desired pressure is indicated on the manometer while taking care not to exceed the desired pressure.
2. To decrease pressure, release force against the plunger thumb ring.
3. To create negative pressure or aspirate, pull back on the plunger thumb ring.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

US Patent No. 5,713,242, 6,796,959B2, and 6,938,319B2

Atrion and QL® are Registered Trademarks of Atrion Medical Products, Inc.

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