QL[®] Inflation Device

QL[®] Inflation Device

Dispositif de gonflage QL[®]

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Dispositivo de inflado QL[®]

Dispositivo di gonfiaggio QL[®]

QL[®] Vulhulpmiddel

QL[®] inflationsanordning

QL[®] fyllningsenhet

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Atrion Medical Products, Inc.
1426 Curt Francis Road
Arab, Alabama 35016 USA

Symbol Legend, French (European), German, Spanish (European), Italian, Dutch, Danish, Swedish, Greek, Portuguese (European), Czech, Hungarian, Polish, Slovak, Norwegian, Finnish, Bulgarian, Romanian, Estonian, Lithuanian, Latvian, Turkish, Croatian, Russian, Japanese, Simplified Chinese

	Contents Contenu Inhalt Contenido Contenuto Inhoud Innhold Innehåll Περιεχόμενα Conteúdo Obsah Tartalom Spis treści	Obsah Innhold Sisältö Содержание Cuprins Sisu Turinys Saturis İçerik Sadržaj Содержание Tartalom 目次 内容物		Manufacturer Produsent Hersteller Fabricante Produttore Fabrikant Producent Tillverkare Κατασκευαστής Fabrico Výrobce Gyártó Producent	Výrobca Produsent Valmistaja Производитель Producător Tootja Gamintojas Ražotājs Üretici Proizvodáč Изготовитель 製造業者 制造商		Lot Number Número de lot Chargennummer Número de lote Numero di lotto Partijnummer Lotnummer Lotnummer Αριθμός παρτίδας Número de lote Číslo šarže Tételszám Numer partii	Číslo šarže Partinummer Eränumero Partiden номер Număr lot Partii number Lotnummer Partijas numurs Parti Numarası Broj serije Номер партии ロット番号 批号		Unique Device Identifier Identifiant unique du dispositif Einmalige Produktkennung Número de dispositivo único Numero univoco del dispositivo Unieke hulpmiddelidentificatie Unikt udstyrsnummer Unikt enhetsnummer Αποκλειστικός αριθμός τεχνολογικού Número único do dispositivo Jedinečné číslo prostředku Egyedi eszközazonosító (UDI) Unikalny numer urządzenia	Jedinečné identifikačné číslo pomůcky Unikt utstyrsnummer Yksilöllinen laitenumero Уникален номер на изделието Număr unic al dispozitivului Kordumatu seadme number Unikalusis priemonės numeris Ierices unikālais numurs Benzersiz Cihaz Tanımlayıcı Jedinstveni broj uređaja Уникальный номер изделия 固有のデバイス番号 器械唯一标识
	Non-pyrogenic Apyrogenė Nicht pyrogen Apirógeno Apirogeno Niet-pyrogeen Ikke-pyrogen Icke-pyrogen Μη πυρετογόνο Apirogénico Nepyrógenni Nem pirogén Niepirogenny	Nepyrógénne Ikke-pyrogen Ei-pyrogeninen Непирогенно Apirogen Mittepyrogeen Nepirogeniškas Nepirógens Pirojenik değıldir Nepirogeno Nepyrógenni 非発熱性 无热原		Date of manufacture Date de fabrication Herstellungsdatum Fecha de fabricación Data di fabbricazione Fabricagedatum Fremstillingsdato Tillverkningsdatum Ημερομηνία κατασκευής Data de fabrico Дата изготовления Gyártás dátuma Data produkci	Dátum výroby Produksjonsdato Valmistuspäivämäärä Дата на производство Data fabricației Tootmiskuupäev Pagaminimo data Izgatavošanas datums Üretim tarihi Datum proizvodnje Дата изготовления 製造日 制造日期		Model Number Référence de modèle Modellnummer Número de modelo Numero di modello Modelnummer Modelnummer Modellnummer Αριθμός μοντέλου Número do modelo Číslo modelu Tipusszám Numer modelu	Číslo modelu Modellnummer Modellnummer Номер на модел Număr de model Mudeli number Modelo numeris Modela numurs Model Numarası Broj modela Номер модели モデル番号 型号		Consult instructions for use Identifiant unique du dispositif Gebrauchsanweisung beachten Consultar las instrucciones de uso Consultare le istruzioni per l'uso Raadpleeg de gebruiksaanwijzing Se brugsanvisningen Se bruksanvisningen Συμβουλευτείτε τις οδηγίες χρήσης Consultar as instruções de utilização Přečtěte si návod k použití Olvasa el a használati utasítást Zapoznać się z instrukcją użytkowania	Přečítajte si návod na použitie Se bruksanvisningen Lue käyttöohjeet Викйте инструкциите за употреба Consultati instrucțiunile de utilizare Žr. naudojimo instrukcijas Skatiet lietošanas pamācību Kullanma talimatına başvurun Pogledajte upute za uporabu См. инструкцию по применению 使用説明書を参照してください 查阅使用说明
	Do not reuse Ne pas réutiliser Nicht wiederverwenden No reutilizar Non riutilizzare Niet opnieuw gebruiken Mā ikke genbruges Får inte återanvändas Μην επαναχρησιμοποιείτε Nāo reutilizar Nepoužívejte opakovaně Ne használja fel újra Nie używać ponownie	Nepoužívejte opakovane Mā ikke gjenbrukes Ei saa käyttää uudelleen Да не се използва повторно A nu se reutiliza Mitte korduvalt kasutada Nenaudoti pakartotinai Nelietot atkārtoti Tekrar kullanmayın Nāo reutilizar Nemojte ponovno upotrebljavati Не использовать повторно 再使用しないでください 请勿重复使用		Do not re-sterilize Ne pas restériliser Nicht reesterilisieren No reesterilizar Non ristertilizzare Niet opnieuw steriliseren Mā ikke resteriliseres Får inte omsteriliseras Μην επαναποστειρώσετε Nāo reesterilizar Neresterilizirati Ne sterilizálja újra Nie sterylizować ponownie	Nesterilizujte opakovane Mā ikke resteriliseres Ei saa steriloida uudelleen Да не се стерилизира повторно A nu se resteriliza Mitte uuesti steriliseerida Nesterilizuoti pakartotinai Nesterilizēt atkārtoti Tekrar sterilize etmeyin Nemojte ponovno sterilizirati Не стерилизовать повторно 再滅菌しないでください 请勿重新灭菌		Medical Device Dispositif médical Medizinprodukt Producto sanitario Dispositivo medico Medisch hulpmiddel Medicinsk udstyr Medicinteknisk produkt Ιατροτεχνολογικό προϊόν Dispositivo médico Zdravotnický prostředek Orvostechnikai eszköz Wyrob medyczny	Zdravotnícka pomôcka Medisinsk utstyr Lääkinnällinen laite Μεδικινησκό изделие Dispozitiv medical Meditsiiseade Medicinos priemonė Medicinskā ierice Tibbi Cihaz Medicinski proizvod Медицинское изделие 醫療機器 医疗器械		Quantity Quantité Menge Cantidad Quantità Aantal Antal Antal Ποσότητα Quantidade Množství Mennyiség Ilość	Množstvo Antall Määrä Количество Capitate Kogus Kiekis Daudzums Miktar Količina Количество 数量 數量

Symbol Legend, French (European), German, Spanish (European), Italian, Dutch, Danish, Swedish, Greek, Portuguese (European), Czech, Hungarian, Polish, Slovak, Norwegian, Finnish, Bulgarian, Romanian, Estonian, Lithuanian, Latvian, Turkish, Croatian, Russian, Japanese, Simplified Chinese

Distributor
Distributeur
Vertriebshändler
Distribuidor
Distributore
Distributeur
Distributor
Distributör
Διανομέας
Distribuidor
Distributor
Forgalmazó
Dystrybutor

Distribütör
Distributeur
Jakelija
Дистрибутор
Distributor
Edasimüüja
Platintojas
Izplatītājs
Distribütör
Distributer
Дистрибутор
販売代理店
经销商


Importer
Importateur
Importeur
Importador
Importatore
Importeur
Importør
Importör
Εισαγωγέας
Importador
Dovozce
Importör
Importer

Dovozca
Importor
Maahantuoja
Вносител
Importator
Importija
Importuotojas
Importētājs
Íthalatçı
Uvoznik
Импортер
輸入者
进口商

Use-by date
Date limite d'utilisation
Verwendbar bis
Fecha de caducidad
Data di scadenza
Uiterste gebruiksdatum
Utlöbsdato
Utgångsdatum
Ημερομηνία λήξης
Prazo de validade
Datum použitelnosti
Lejárató dátum
Data ważności

Dátum použitelnosti
Utlöbsdato
Viimeinen käyttöpäivämäärä
Дата на изтичане на срока на годност
Data de expirare
Aegumiskuupäev
Galiojimo pabaigos data
Derīguma termiņš
Son kullanma tarihi
Rok valjanosti
Срок годности
使用期限
有效期

Tubing length
Longueur de la tubulure
Schlauchlänge
Lengitud del tubo
Lunghezza del tubo
Slanglengte
Slangelængde
Slanglängd
Μήκος σωλήνωσης
Comprimento da tubagem
Délka hadičky
Cső hossza
Długość rurki

Dĺžka hadičky
Slangelengde
Lektun pituus
Дължина на тръбата
Lungimea tubulaturii
Vooliku pikkus
Vamzdelio ilgis
Skyddsförpackning med en steril barriär på utsidan
Προστατευτική συσκευασία μονού στείρου φραγμού εξωτερικά
Embalagem protetora de barreira esteril única no exterior
Jedno sterilní bariérové ochranné balení vně
Egyszeres steril védőzáras védőcsomagolás kívül
Opakowanie ochronne z pojedynczą barierą sterylą od zewnątrz

Sterilized using ethylene oxide
Stérilisé à l'oxyde d'éthylène
Sterilisation mit Ethylenoxid
Esterilizado con óxido de etileno
Sterilizato con ossido di etilene
Gesteriliseerd met ethyleenoxide
Steriliseret med ethylenoxid
Steriliserad med etylenoxid
Αποστειρωμένο με οξείδιο του αιθyleνίου
Esterilizado por óxido de etileno
Sterilizováno ethylenoxidem
Etilén-oxidál sterilizálva
Steryizowano tlenkiem etylenu

Sterilizované etylénoxidom
Steriliseret med etylenoksid
Steriloitu eteenoksidilla
Стерилизовано с этиленов оксид
Sterilizat cu oxid de etilenă
Steriliseeritud etüleenoksiidiga
Sterilizuota etileno oksidu
Sterilizēts, izmantojot etilēnoksidu
Etilen oksit kullanılarak sterilize edilmiştir
Sterilizirano etilen-oksidom
Стерилизовано оксидом этилена
エチレンオキシド滅菌済み
已使用环氧乙烷灭菌

Single sterile barrier protective packaging outside
Emballage protecteur à barrière stérile unique à l'extérieur
Einzelne Sterilbarriere mit äußerer Schutzverpackung
Embalaje protector de barrera esteril única exterior
Confezione protettiva singola barriera sterile esterna
Toepassing van een enkelvoudige steriele barrière binnen een beschermende verpakking
Beskyttende emballage med enkelt steril barriere udvendigt
Skyddsförpackning med en steril barriär på utsidan
Προστατευτική συσκευασία μονού στείρου φραγμού εξωτερικά
Embalagem protetora de barreira esteril única no exterior
Jedno sterilní bariérové ochranné balení vně
Egyszeres steril védőzáras védőcsomagolás kívül
Opakowanie ochronne z pojedynczą barierą sterylą od zewnątrz

Ochranný obal s jedinou sterilnou bariérou na vonkajšej strane
Beskyttelsemballasje med enkel steril barriere på utsiden
Yksittäinen steriiliesteen suojaopakkaus ulkopuolella
Единична стерилна защитна опаковка отвън
Ambalaj de protecție cu barieră sterilă unică în exterior
Väljaspool ühekordne steriilse barjääriga kaitsepakend
Vieno sterilumo barjero apsauginė pakuotė išorėje
Vienas sterilās barjeras aizsargiepakojums ārpusē
Diş kısmıda tek steril bariyer koruyucu ambalaj
Vanjska zaštitna ambalaža s jednostrukom sterilnom barijerom
Одноразный стерильный барьер, защитная упаковка снаружи
外側が保護包装されているシングル無菌バリヤ
单层无菌屏障保护性外包装

Do not use if package is damaged and consult instructions for use
Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi
Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten
No utilizar si el envase está dañado y consultar las instrucciones de uso
No utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
Niet gebruiken als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing
Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen
Använd inte om förpackningen är skadad och läs bruksanvisningen
Μη χρησιμοποιείτε, εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης
Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
Nepoužívejte, pokud je obal poškozený, a prostudujte si návod k použití
Ne használni, ha a csomagolás sérült, és olvassa el a használati utasítást
Nie używać, jeśli opakowanie jest uszkodzone i zapoznać się z instrukcją użytkowania

Nepoužívejte, ak je obal poškodený a prečítajte si návod na použitie
Må ikke brukes hvis pakningen er skadet, og se bruksanvisningen
Ei saa käyttää, jos pakkaus on vaurioitunut, ja katso käyttöohjeet
Не използвайте, ако опаковката е повредена, и направете справка в инструкциите за употреба
A nu se utiliza dacă ambalajul este deteriorat și consultați instrucțiunile de utilizare
Ärge kasutage, kui pakend on kahjustatud, ja vt kasutusjuhendit
Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas pamācību
Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun
Nemojte koristiti ako je pakiranje oštećeno i pročitajte upute za uporabu
Не использовать, если упаковка повреждена, и обратиться к инструкции по применению
包装が破損している場合は使用せず、取扱説明書を参照してください。
包装如有破损，请勿使用，并查阅使用说明

Single sterile barrier system
Système de barrière stérile unique
Einzelnes Sterilbarriersystem
Sistema de barriera esteril única
Sistema a barriera sterile singola
Toepassing van een enkelvoudige steriele barrière
Sterilt enkeltbarriersystem
Enkelt steril bariérsystem
Σύνστημα μονού στείρου φραγμού
Sistema de barreira esteril única
Systém jednoduché sterilní bariéry
Egyszeres steril védőzáras rendszer
System pojedynczej bariery sterylnej

Systém jedinej sterilnej bariéry
System med enkel steril barriere
Yksittäinen steriiliestearjestelmä
Система с единична стерилна барьера
Sistem cu barieră sterilă unică
Ühekordne steriilse barjääri süsteem
Vieno sterilumo barjero sistema
Vienas sterilās barjeras sistēma
Tek steril bariyer sistemi
Sustav jednostruke sterilne barijere
Одноразная стерильная барьерная система
シングル無菌バリアシステム
单层无菌屏障系统


Authorized representative in the European Union
Mandataire établi dans l'Union européenne
Bevollmächtigter in der Europäischen Union
Representante autorizado en la Unión Europea
Rappresentante autorizzato nell'Unione Europea
Gemachtigde in de Europese Unie
Autoriseret representant i EU
Auktoriserad representant i Europeiska unionen
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση
Representante autorizado na União Europeia
Zplnomocnený zástupce v Európskej únii
Meghatalmazott képviselő az Európai unióban
Autoryzowany przedstawiciel w Unii Europejskiej

Splnomocnený zástupca v Európskej únii
Autoriseret representant i Den europeiske union
Valtuutettu edustaja Euroopan unionissa
Упълномощен представител в Европейския съюз
Reprezentant autorizat în Uniunea Europeană
Volitatud esindaja Euroopa Liidus
Igalotasis atstovas Europos Sąjungoje
Pilnvarotais pārstāvis Eiropas Savienībā
Avrupa Birliği'ndeki yetkili temsilci
Ovlašteni zastupnik u Europskoj uniji
Уполномоченный представитель в Европейском союзе
歐洲連合の正式な代表者
歐盟授權代表

Store in a cool, dry place
Conserver dans un endroit frais et sec
Kühl und trocken lagern
Almacenar en un lugar fresco y seco
Conservare in luogo fresco e asciutto
Op een koele, droge plaats bewaren
Opbevares køligt og tørt
Förvaras på en sval och torr plats
Φυλάσσετε σε δροσερό, ξηρό χώρο
Armazene num local fresco e seco
Skladujte na chladném, suchém místě
Hüvös, száraz helyen tárolandó
Przechowywać w chłodnym i suchym miejscu

Skladujte na chladnom a suchom mieste
Oppbevares på et kjølig, tørt sted
Säilytetävä viileässä ja kuivassa paikassa
Да се съхранява на хладно, сухо място
A se depozită într-un loc răcoros și uscat
Hoida jahedas ja kuivas kohas
Laikyti vėsioje, sausoje vietoje
Uzglabāt vēsā, sausā vietā
Serin, kuru bir yerde saklayın
Čuvajte na hladnom i suhom mjestu
Хранить в сухом прохладном месте
涼しく乾燥した場所に保管してください
存放在阴凉干燥处

QL® Inflation Device
Dispositif de gonflage QL®
QL® Aufdehnungsvorrichtung
Dispositivo de inflado QL®
Dispositivo di gonfiaggio QL®
QL® Vulhulpmiddel
QL® inflationsanordning
QL® fyllningsenhet
Συσκευή διόγκωσης QL®
Dispositivo de insuflação QL®
Plnicí zařízení QL®
QL® feltöltőeszköz
Urządzenie do napiewiania QL®

Plniaca pomôcka QL®
QL®-fyllingsenhet
QL®-täyttölaitte
Устройство за раздуване QL®
Dispozitiv de umflare QL®
QL® täitmiseseade
QL® isplétimo prietaisas
QL® uzpildes ierice
QL® şişirme cihazı
Uredaj za napuhavanje QL®
Устройство для накачивания баллона QL®
QL® インフレーションデバイス
QL® 充盈装置

Stopcock
Robinet
Absperrhahn
Llave de paso
Rubinetto
Afsluitkraan
Stopfane
Kran
Στρόγγυς
Torneira de passagem
Uzavírací kohout
Zárócsap
Zawór

Uzatvárací kohútik
Stoppekran
Sulkuhana
Спирателен кран
Robinet
Sulgurkraan
Čiaupas
Noslēgkrāns
Ūglī valf
Sigurnosni ventil
Стопорный кран
止栓
旋塞閥

Unlock/Lock
Déverrouillage/verrouillage
Entriegeln/Verriegeln
Desbloquer/bloquer
Sblocca/blocca
Ontgrendelen/Vergrendelen
Läs op/läs
Läs upp/läs
Αποσφάλιση/ασφάλιση
Desbloquer/bloquer
Odemknout/zamknout
Kioladás/zárolás
Odblokowanie/zablokowanie

Odomknút/uzamknút
Uläst/läst
Avaa/lukitse
Отключване/заклучване
Debloccare/bloccare
Lukust avamine / lukustamine
Atlaisvinti / užfiksuoti
Atbloķēšana/bloķēšana
Kilit Aç/Kilitle
Otključavanje/zaključavanje
Разблокировать/заблокировать
ロック解除／ロック
解锁/锁定

Rx only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

QL® Inflation Device - Instructions for Use

The product is intended for use by physicians trained and experienced in interventional techniques.

Device Description:

The QL® Inflation Device is a disposable single use catheter balloon inflator with a lock lever design that controls the piston plunger, a pressure gauge, and a connecting tube with a male luer lock connector. The pressure gauge measures pressures ranging from vacuum to gauge capacity; the gauge is marked in 1 atm increments. The gauge may also have an inner scale of comparable PSI measurements. The accuracy of the pressure gauge has been determined to be per the following table:

Maximum Gauge Atmosphere (atm)	Gauge Accuracy
Up to 30	±1 atm over the range
31 to 55	±1.65 atm over the range

Intended Use:

The QL® Inflation Device is intended to inflate or deflate balloon interventional devices.

Indications:

The QL® Inflation Device is recommended for use while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon.

Contraindications:

There are no known contraindications

Warnings:

- Use only liquid inflation media. Do not inflate with air.
- Before use, confirm that the gauge responds to pressure freely.
- Always follow the manufacturer's directions accompanying the balloon catheter for instructions for use, maximum balloon inflation pressure, precautions, and warnings for that device.
- The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- EU MDR 2017/745 Annex I 10.4 does not apply to this type of device per Annex I 10.4.1 device descriptions.

Potential Adverse Events:

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

Precaution:

- Before use, inspect the device to verify that no damage has occurred during shipping and handling.
- Before use, ensure the inflation device and tubing is completely free of air.

Directions for Use:

Preparation:

Make all aspiration and injection maneuvers with the lock lever unlocked. Unlock the piston plunger by pushing the lock lever to the unlocked position. In this position, you can freely pull the piston plunger back for aspiration, or push it forward for injection. To lock the piston plunger in position, place the lock lever in the straight up position.

1. Prepare a solution of contrast media and/or normal saline in a small sterile bowl. Check balloon interventional device and contrast media instructions for specific contrast mixture recommendations.
2. Orient the tubing end downward into the contrast media solution.
3. Withdraw unlocked piston plunger to aspirate enough solution to fill the syringe. Attach stopcock and ensure it is in open position.
4. Hold the device upright to purge the air from the syringe, connecting tube, and stopcock. Tap the syringe lightly, if necessary, advance the piston plunger to push all air bubbles out through the tubing and stopcock.

5. Inspect the inflation device, tubing, and stopcock to ensure that the device has been completely purged of air bubbles. If not, repeat steps 2-4 (without removing the stopcock) until no air bubbles are present in the fluid path.
6. Adjust the inflation device volume to the desired amount. If more contrast media solution and/or saline is needed, submerge the tubing end downward into the basin of solution and adjust plunger piston position, lock the lever and close stopcock.
7. Generate pressure by turning the piston plunger clockwise and verify gauge needle movement to confirm the gauge's pressure mechanism and/or movement suffered no damage in transit/handling.

Attaching the inflation device to the balloon dilation catheter:

1. Prepare and test the balloon interventional device according to the directions for use.
2. If a separate syringe was used to prepare the balloon interventional device, remove it. Create a fluid-fluid connection between the stopcock and balloon catheter placing a drop of contrast solution from a syringe into each hub.
3. Hand-tighten the hubs securely.

Operating inflation device:

1. Fill device with the recommended amount of solution, then place the lock lever in the locked straight up position to allow fine volume adjustments or maintain the piston at the desired volume setting.
2. To rapidly inject fluid into the balloon, release the lock lever and manually push the plunger forward until you reach the balloon manufacturer's recommended volume. To generate pressure in the balloon, lock the lever and rotate piston plunger clockwise
3. For rapid deflation, with the lock lever in its unlocked position, pull back the piston plunger manually.

How Supplied:

The QL® Inflation Device is supplied in a tray with a Tyvek lid and sterilized with ethylene oxide.

Device Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

Hazardous Substances:

This medical device does not contain hazardous substances.

The following US Patent No. covers the Atrion Inflation Device: 10,390,868

▲**Atrion** and QL® are Registered Trademarks of Atrion Medical Products, Inc.

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THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK