Endoscopic Applicator

INSTRUCTIONS FOR USE

DESCRIPTION: The Endoscopic Applicator is a sterile, single-use applicator device. The Endoscopic Applicator consists of two components: (1) a non-reflective composite cannula, and (2) a composite stylet.

INDICATIONS: The Endoscopic Applicator is intended for use in delivering hemostatic agents to bleeding surgical sites through a 5 mm or larger trocar.

HOW SUPPLIED: Each Endoscopic Applicator is supplied sterile in a sealed plastic container. One (1) box contains 6 Endoscopic Applicators.

WARNINGS/PRECAUTIONS:

- Read the hemostatic agent instructions for use prior to using this device.
- Do not use if damaged.
- Do not use air as a propellant to extrude hemostatic agents instead of the stylet.
- Do not use Endoscopic Applicator to manipulate or retract organs or tissue.
- The Endoscopic Applicator should only be used by persons adequately trained in Minimal Invasive procedures.
- Do not reuse or resterilize. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

DIRECTIONS FOR USE:

1. Single use only – see warnings/precautions.
2. Read the hemostatic agent instructions for use prior to using this device.
3. Transfer both cannula and stylet (obturator) to the sterile field.
4. Attach the syringe with the hemostatic agent to the cannula, ensuring the luer connector is locked.
5. Load the Endoscopic Applicator with the desired amount of hemostatic agent by depressing the plunger on the hemostatic agent syringe. The cannula can be filled with a maximum of 5ml of hemostatic agent.
6. Introduce the Endoscopic Applicator cannula into the trocar port. **CAUTION:** Grasp the cannula while advancing it through the trocar to avoid depressing the syringe plunger.
7. Carefully position the distal end of the Endoscopic Applicator at the site where the hemostatic agent is to be delivered.
8. Detach the hemostatic agent syringe from the luer lock connector of the Endoscopic Applicator.
9. Use the stylet to dispense the hemostatic agent to the bleeding site in accordance with the instructions for use of the hemostatic agent.
10. When sufficient hemostatic agent has been delivered to the bleeding site, carefully remove the Endoscopic Applicator from the trocar port, remove the stylet from the cannula and flush any remaining hemostatic agent from the cannula.
11. If additional hemostatic agent is required during the same surgical procedure, repeat the steps above as necessary.

Single-use, one patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of patient.

United States only: 1-800-423-2090
Other Countries: Contact your local Baxter representative.

USA: Rx only

Baxter is a trademark of Baxter International Inc.

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Baxter er et varemærke, der tilhører Baxter International Inc.