Malleable Tip Endoscopic Applicator

INSTRUCTIONS FOR USE

DESCRIPTION: The Malleable Tip Endoscopic Applicator is a sterile, single-use device that consists of two components: (1) a cannula (composite tube with a malleable tip), and (2) a composite stylet.

INDICATIONS: The Malleable Tip Endoscopic Applicator is intended for use in delivering a hemostatic agent to bleeding sites through a 5 mm or larger diameter trocar.

HOW SUPPLIED: Each Malleable Tip Endoscopic Applicator is supplied sterile in a sealed plastic container. One (1) box contains six (6) Malleable Tip Endoscopic Applicators.

WARNINGS/PRECAUTIONS:

- Read the hemostatic agent instructions for use prior to using this device.
- Do not use if damaged.
- Do not use air as a propellant to extrude hemostatic agents instead of the stylet.
- Do not use Malleable Tip Endoscopic Applicator to manipulate or retract organs or tissue.
- Do not bend the malleable tip against organ or tissue.
- The Malleable Tip Endoscopic Applicator should only be used by persons adequately trained in Minimally Invasive surgical procedures.
- Do not reuse or resterilize. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

DIRECTIONS FOR USE:

1. Single use only – see warnings/precautions.
2. Read the hemostatic agent instructions for use prior to using this device.
3. Transfer both cannula and sylet (obturator) to the sterile field.
4. Attach the syringe with the hemostatic agent to the cannula, ensuring the Luer connector is locked.
5. Load the Malleable Tip Endoscopic Applicator with the desired amount of hemostatic agent by depressing the plunger on the hemostatic agent syringe. The cannula can be filled with a maximum of 5.8mL of hemostatic agent.
6. Ensure the malleable tip is straight before introducing the Malleable Tip Endoscopic Applicator cannula into the trocar port. **NOTE:** Grasp the cannula while advancing it through the trocar to avoid depressing the syringe plunger.
7. Carefully approximate the distal end of the Malleable Tip Endoscopic Applicator to the treatment site. Use appropriate laparoscopic instruments to bend the malleable tip.
8. Detach the hemostatic agent syringe from the luer lock connector of the Malleable Tip Endoscopic Applicator.
9. Use the stylet to dispense the hemostatic agent to the bleeding site in accordance with the instructions for use of the hemostatic agent.
10. After application of the hemostat, straighten the malleable tip with appropriate laparoscopic instruments, carefully remove the Malleable Tip Endoscopic Applicator from the trocar port, remove the stylet from the cannula, and flush any remaining hemostatic agent from the cannula.
11. If additional hemostatic agent is required during the same surgical procedure, repeat the steps above as necessary.
12. Properly dispose of the Malleable Tip Endoscopic Applicator as Biohazard waste.

Single-use, one patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of patient.

USA: Rx only

Baxter is a trademark of Baxter International Inc.

United States only: 1-800-423-2090

Other Countries: Contact your local Baxter representative.

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Baxter je ochranná známka společnosti Baxter International Inc.