
Baxter

DuploTip

DUAL LUMEN CANNULA

Instructions for Use
Mode d'emploi
Gebruiksaanleitung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Указания за употреба
Pokyny k použití
Brugsanvisning
Οδηγίες Χρήσης
Käyttöohjeet
Uputstva za uporabu
Instrukcja obsługi
Bruksanvisning
Kullanım Talimatları

DIST

Baxter Healthcare Co.
Deerfield, IL 60015, USA
1-888-229-0001

Manufactured in Mexico for:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977

CE 0297

EC

REP

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

PUSH HERE
TO RELEASE

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Dual Lumen Cannula

Before using, read the following information:

DESCRIPTION / INDICATIONS FOR USE:

Sterile, single-use device for the application of TISSEEL Two-Component Fibrin Sealant

CONTRAINDICATIONS:

None

INSTRUCTIONS FOR USE:

1. Remove the Applicator from its package using sterile technique. Remove the protective tube.
2. Push the male syringe set Luers firmly into the female Luer cones on the Applicator as shown in Figure 1.1.
3. Push down on the snap lock to securely fasten the Applicator to the syringe Luers as shown in Figure 1.2.
4. Tips (metal only) may be angled to improve access or visibility of surgical site.
5. Applicator is ready for use.
6. To release the Applicator from the syringe Luers, press the release button on the Snap Lock as shown in Figure 2.
7. Dispose of as biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient.

IMPORTANT NOTE:

This system consists of two parts, a Tisseel syringe system and a dual cannula Applicator. Tisseel syringe systems and applicators must be ordered separately. Make sure that you have both parts before proceeding.

Baxter, DuploTip and Tisseel are trademarks of Baxter International Inc., or its subsidiaries.

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Baxter, DuploTip ve Tisseel; Baxter International Inc. veya bağlı şirketlerinin ticari markalarıdır.
Baxter, DuploTip och Tisseel är varumärken som tillhör Baxter International Inc. eller dess dotterbolag.