www.nordsonmedical.com/efiu

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Istruções de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização
Kullanma Talimatı

MANUFACTURED FOR:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977
800-624-5662
info@nordsonmedical.com
www.nordsonmedical.com



CE 0297

EC REP

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FibriJet® APPLICATOR ASSEMBLY

Before using, read the following information:

****IMPORTANT NOTE****
This device is intended for use with an Applicator Tip which is not included. Applicator Tips must be ordered separately. Make sure you have both parts before proceeding.

DESCRIPTION / INDICATIONS FOR USE:

Sterile, single-use instrument for the application of two liquids.

CONTRAINDICATIONS:

None known.

WARNINGS / PRECAUTIONS:

Device is intended for use by Medical Professionals.

If sterile packaging is damaged or unintentionally opened prior to use, do not use the device and dispose per the Disposal instructions below. Device is intended for single use only. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient. There is a low probability that a delay in treatment, increased procedural time, mild tissue damage or local immunogenic response could occur if the device fails to function as intended.

INSTRUCTIONS FOR USE:

- REMOVE the Applicator Tip (not included) from its package and set aside.
- REMOVE the Applicator Assembly from its sterile package using sterile technique.
- FILL the two Applicator syringes with desired volumes of liquids [3ml 0-1±.07, 1.0-3.0±.05; 5ml 0-3±.18, 3.0-5.0±.05; 10ml 0-4.0±.18, 4-10±.48].
- (A) ASSEMBLE the syringes to FibriJet Applicator Tip by firmly twisting syringes onto luer-lock hubs of Tip.

- (B) ASSEMBLE the holder and plunger clip to syringes as shown.
- (C) Applicator is ready for use.
- The Applicator may be refilled for subsequent uses **in the same procedure only.**
To Refill Applicator follow procedure below:
 - Remove the plunger link.
 - Remove syringes from holder.
 - Remove the Applicator Tip from the syringes by twisting both syringes counter clockwise.
 - Refill the syringes.
 - Reinstall the Applicator Tip as described above.**IMPORTANT: To avoid clogging applicator, always connect larger syringe to port with the round indicator.**

DISPOSAL:

Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

STORAGE CONDITIONS:

Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity.


Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member state in which the user/patient is established.

United States only: 1-800-624-5662

Other Countries: Contact your local representative and/or distributor

INTERNATIONAL LANGUAGE SERVICES
601 CARLSON PKWY #1050
MINNETONKA, MN 55305
800-225-8964

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