
FibriJet®
BLENDING CONNECTOR

www.nordsonmedical.com/eifu

Instructions for Use
Mode d'emploi
Gebruiksaanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização

MANUFACTURED FOR:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977
800-624-5662



CE 0297

EC REP

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FibriJet®

BLENDING CONNECTOR

Before using, read the following information:

******IMPORTANT NOTE******
This device is intended for use with an
Applicator Assembly which is not included.
Applicator Assemblies must be ordered
separately. Make sure you have both parts
before proceeding.

DESCRIPTION / INDICATIONS FOR USE:
Sterile, single-use instrument for the application of two liquids.

CONTRAINDICATIONS:
None known.

WARNINGS / PRECAUTIONS:
If sterile packaging is damaged or unintentionally opened prior to use, do not use the device and dispose per the Disposal instructions below.

Device is intended for single use only. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient.

There is a low probability that a delay in treatment, increased procedural time, mild tissue damage or local immunogenic response could occur if the device fails to function as intended.

- INSTRUCTIONS FOR USE:**
- REMOVE Applicator Tip from its sterile package using sterile technique.
 - REMOVE the Applicator Assembly from its sterile package.
 - See Applicator Assembly Instructions (**FIG. 1**).
 - ASSEMBLE included Spray Tip or Cannula Tip to Blending Connector by twisting firmly clockwise onto luer-lock connection (**FIG. 2**).
 - COMPLETE ASSEMBLY according to Applicator Assembly Instructions.
 - Applicator Tip is ready for use.

NOTE: To avoid clogging, always put syringe containing Fibrinogen or Plasma on side with round indicator.

DISPOSAL:
Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

STORAGE CONDITIONS:
Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member state in which the user/patient is established.
United States only: 1-800-624-5662
Other Countries: Contact your local representative and/ or distributor

INTERNATIONAL LANGUAGE SERVICES
601 CARLSON PKWY #1050
MINNETONKA, MN 55305
800-225-8964

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