
FibriJet®
Single Spray Tip

www.nordsonmedical.com/EIFU

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização

Manufactured For:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977


CE 0297

EC REP

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EN

FibriJet®
Single Spray Tip

Before using, read the following information:

****IMPORTANT NOTE****
This device is intended for use with a syringe or compatible applicator with a standard luer lock which is not included and must be ordered separately. Make sure you have both parts before proceeding.

DESCRIPTION/ INDICATIONS FOR USE
Sterile, single-use spray tip for the application of liquid(s).
CONTRAINDICATIONS:
None known

WARNINGS / PRECAUTIONS:
Ensure spray tip is secured to prevent accidental tip detachment.
Device is intended for use by Medical Professionals.
If sterile packaging is damaged or unintentionally opened prior to use, do not use the device and dispose per the Disposal instructions below.
Device is intended for single use only. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient.
There is a low probability that a delay in treatment, increased procedural time, mild tissue damage or local immunogenic response could occur if the device fails to function as intended.

INSTRUCTIONS FOR USE:

- REMOVE the spray tip from its sterile package using sterile technique.
- (A) Twist spray tip onto a standard luer lock (i.e., syringe)
-Ensure spray tip is secured to prevent accidental tip detachment.
- (B) Applicator tip is ready to use.

DISPOSAL
Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

STORAGE CONDITIONS:
Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member state in which the user/patient is established.
United States only: 1-800-624-5662
Other Countries: Contact your local representative and/or distributor

INTERNATIONAL LANGUAGE SERVICES
601 Carlson Pkwy #1050
Minnetonka, MN 55305
800-225-8964

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