Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Modo de empleo
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning

Manufactured For:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977
800-624-5662



CE 0297

EC REP

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EN

FibriJet® Mixing Applicator Tips

Before using, read the following information:

****IMPORTANT NOTE****
This device is intended for use with an Applicator Assembly which is not included. Applicator Assemblies must be ordered separately. Make sure you have both parts before proceeding.

DESCRIPTION / INDICATIONS FOR USE:
Sterile, single-use instrument for the application of two liquids.

CONTRAINDICATIONS:
None known.

WARNINGS / PRECAUTIONS:
Device is intended for use by Medical Professionals.
If sterile packaging is damaged or unintentionally opened prior to use, do not use the device and dispose per the Disposal instructions below.
Device is intended for single use only. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient.
There is a low probability that a delay in treatment, increased procedural time, mild tissue damage or local immunogenic response could occur if the device fails to function as intended.


- INSTRUCTIONS FOR USE:**
- REMOVE Applicator Tip from its sterile package using sterile technique.
 - REMOVE the Applicator Assembly from its sterile package using sterile technique.
 - Twist Applicator Assembly syringes firmly onto luer-lock hubs of Mixing Tip. (See Applicator Assembly Instructions) **(FIG 1)**.
 - ASSEMBLE included Spray Tip to Mixing Tip by twisting the rotating collar clockwise onto luer-lock connection **(FIG 2.1)**. **Applicators for high viscosity fluids do not contain the Spray Tip (FIG 2.2)**.
 - COMPLETE ASSEMBLY according to Applicator Assembly Instructions.
 - Applicator Tip is ready for use.

DISPOSAL:
Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

STORAGE CONDITIONS:
Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member state in which the user/patient is established.
United States only: 1-800-624-5662
Other Countries: Contact your local representative and/or distributor

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