Instructions for Use

Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização

56 EAST BELL DRIVE
P.O. BOX 587
WARSAW, IN 46581 USA

Manufactured for:

Micromedicals, Inc.
1270 Eagan Industrial Road
St. Paul, MN 55121-1385 USA
651-452-1977

CE 0297

EC REP

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

B

EN
MIXING CONNECTOR WITH SPRAY TIP

Before using, read the following information:

****IMPORTANT NOTE****

This device is intended for use with an applicator assembly which is not included. Applicator assembly must be ordered separately. Make sure you have both parts before proceeding.

DESCRIPTION/INDICATIONS FOR USE:

Sterile, single-use instrument for application of two liquids.

CONTRAINDICATIONS

None known.

WARNINGS / PRECAUTIONS

IMPORTANT: To avoid clogging applicator, always connect larger syringe to port with the round indicator.

Device is intended for use by Medical Professionals.

If sterile packaging is damaged or unintentionally opened prior to use, do not use the device and dispose per the Disposal instructions below.

Device is intended for single use only. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of a patient.

There is a low probability that a delay in treatment, increased procedural time, mild tissue damage or local immunogenic response could occur if the device fails to function as intended.

INSTRUCTIONS FOR USE

1. REMOVE the Mixing Connector and Spray Tip from its sterile package using sterile technique.
2. (Figure A) ASSEMBLE the syringes to the Mixing Connector by firmly twisting the syringes clockwise onto luer-lock hubs of the tip.
3. (Figure B) ASSEMBLE Spray Tip to the tip of the Mixing Connector by twisting the rotating collar clockwise onto the luer-lock connection.

DISPOSAL

Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

STORAGE CONDITIONS:


Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member state in which the user/patient is established.

United States only: 1-800-624-5662

Other Countries: Contact your local representative and/or distributor.

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