# Arthrex® BioXpress

**Graft Delivery Device**

**Angle Tip**

Instructions for Use  
Mode d'emploi  
Bruksanvisning  
Gebruiksaanwijzing  
Instruções de utilização  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Kullanma Talimatı  
Brugsanvisning  
Instrukcja użytkowania

Manufactured for:

Arthrex, Inc.  
1370 Creekside Blvd.  
Naples, Florida 34108 USA  
(800) 934-4404

Arthrex GmbH  
Erwing-Hielscher-Strasse 9  
81249 München, Germany  
Tel: +49 89 909005-0  
[www.arthrex.de](http://www.arthrex.de)

Made in Mexico

Micromedics, Inc.  
1270 Eagan Industrial Road  
St. Paul MN 55121-1385, USA  
651-452-1977

**CE 0297**

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
## VENTED SYSTEM

| Vented Graft Delivery Device Parts List |                                                   |                           |
|-----------------------------------------|---------------------------------------------------|---------------------------|
| A                                       | Syringe Body                                      |                           |
| B                                       | Plunger Body                                      |                           |
| C                                       | Distal Stylet (contained in Plunger Body)         |                           |
| D                                       | Hydration Tube                                    | Either D or E is present. |
| E                                       | Hydration Adapter                                 |                           |
| F                                       | Vented Stylet Stopper (contained in Plunger Body) |                           |
| G                                       | Plunger Stopper                                   |                           |
| H                                       | Funnel                                            |                           |
| I                                       | Cap Plug (blue)                                   |                           |
| J                                       | Tapered End Cap                                   |                           |
| K                                       | Secondary Thumb Tab                               |                           |
| L                                       | Proximal Stylet (10cm or 15cm Lengths)            |                           |
| M                                       | Angle Tip Cannula (10cm or 15cm Lengths)          |                           |

### 1. INTENDED USE

The Graft Preparation and Delivery Device is intended for the delivery of hydrated allograft, autograft or synthetic bone graft materials to an orthopaedic surgical site. In addition, it is designed to facilitate the premixing of bone graft materials with fluids such as I.V. fluids, blood, plasma concentrate, platelet rich plasma, bone marrow or other specified blood components as deemed necessary by the clinical use requirements.

### 2. WARNINGS/ PRECAUTIONS

There are no specific warnings, precautions, or adverse effects associated with the use of this device.

### 3. INSTRUCTIONS FOR USE

1. Thread Syringe Body (capacity: 15±0.5mL) onto the Tapered End Cap. Thread Hydration Tube or Cap Plug (blue) onto Tapered End Cap.
2. Place Funnel into top of Syringe Body. Place the Tapered End Cap into the retention feature of the tray. Fill Syringe Body with graft material. Remove the Funnel.
3. Push the end of the Distal Stylet into the Plunger Body.
4. Place Plunger Body into Syringe Body. Ensure that Plunger Body has passed the internal retention feature.
5. (A) Device with Hydration Adapter ONLY:  
Rotate to horizontal, then remove Cap Plug and replace with the Hydration Adapter. Hydrate (according to graft material's instructions for use or clinical use requirements) with standard male luer lock syringe (not provided).

#### OR

(B) Device with Hydration Tube ONLY:

Rotate to horizontal, then hydrate (according to graft material's instructions for use or clinical use requirements) using the Hydration Tube with standard male luer lock syringe (not provided).

NOTE: Some material or fluid may pass through the vents in the stopper during hydration.

6. Rotate assembly to horizontal then remove the male luer lock syringe. Remove the Hydration Tube or Hydration Adapter.
7. Thread Cannula onto Tapered End Cap.
8. Depress (push) Plunger Body to dispense material to surgical site.  
NOTE: Some material or fluid may pass through the vents in the stopper during delivery.
9. Thread Proximal Stylet onto Distal Stylet.
10. Depress (push) Stylet through the Cannula to deliver material in the Cannula to the surgical site. Use Secondary Thumb Tab for longer Cannula options.

### 4. DISPOSAL


Dispose of as biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

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