Graft Delivery Device

Instructions for Use
Mode d'emploi
Bruksanvisning
Gebruiksaanwijzing
Instruções de utilização
Gebrauchsanweisung
Istruzioni per l'uso
Instrucciones de uso
Kullanma Talimatı
Brugsanvisning
Instrukcja użytkowania

DIST

Zimmer Biomet Spine, Inc.
10225 Westmoor Dr.
Westminster, CO 80021 USA
+1 800.447.3625

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977

CE 0297

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
CLOSED SYSTEM

Closed Graft Delivery Device Parts List	
A	Syringe Body
B	Plunger Body
C	Distal Stylet (contained in Plunger Body)
D	Closed Stylet Stopper (contained in Plunger Body)
E	Plunger Stopper
F	Funnel
G	Cap Plug (blue)
H	Tapered End Cap
I	Secondary Thumb Tab
J	Proximal Stylet
K	Cannula

1. INTENDED USE

The Graft Preparation and Delivery Device is intended for the delivery of hydrated allograft, autograft or synthetic bone graft material to an orthopedic surgical site. In addition, it is designed to facilitate the premixing of bone graft materials with fluids such as I.V. fluids, blood, plasma concentrate, platelet rich plasma, bone marrow or other specified blood components as deemed necessary by the clinical use requirements.

2. WARNINGS/ PRECAUTIONS

There are no specific warnings, precautions, or adverse effects associated with the use of this device.

3. INSTRUCTIONS FOR USE

1. Thread Syringe Body (capacity: 15±0.5mL) onto the Tapered End Cap. Thread Cap Plug (blue) onto Tapered End Cap.
2. Place Funnel into top of Syringe Body. Place the Tapered End Cap into the retention feature of the tray. Fill Syringe Body with hydrated graft material (according to graft material's instructions for use or clinical use requirements). Remove the Funnel.
3. Push the end of the Distal Stylet into the Plunger Body.
4. Place Plunger Body into Syringe Body. Ensure that Plunger Body has passed the internal retention feature.
5. Rotate to horizontal, then remove Cap Plug.
6. Thread Cannula onto Tapered End Cap.
7. Depress (push) Plunger Body to dispense material to surgical site.
8. Thread Proximal Stylet onto Distal Stylet.
9. Depress (push) Stylet through the Cannula to deliver material in the Cannula to the surgical site. Use Secondary Thumb Tab for longer Cannula options.

4. DISPOSAL


Dispose of as biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

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