
Graft Delivery Device

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização

Manufactured in Mexico for:


Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul, MN 55121- 1385 USA
651-452-1977

CE 0297

EC REP

MDSS GmbH
Schigraben 41
D-30175 Hannover, Germany

Symbol Reference Key / Définition des symboles / Erklärung der Symbole / Clave de referencia a los símbolos / Legenda dei simboli / Verklaring van symbolen / Nøkkel til symbolreferanser / Legenda dos Símbolos		
Graft Delivery Device		Graft Delivery Device Dispositif de pose de greffon Transplantatabgabegerät Dispositivo para suministro de un injerto Dispositivo di erogazione di innesto Hulpmiddel voor graftplaatsing Graftoverføringssystem Dispositivo de Colocação de Enxerto
LATEX	Not made with natural rubber latex Non fabriqué à partir de latex naturel Nicht aus Naturkautschulatelax gefertigt No està fabricado con látex de goma natural Non fabbricato con lattice di gomma naturale Niet vervaardigd met natuurlijke rubberlatex Ikke produsert med lateks av naturgummi Não contém látex de borracha natural	
	55°F - 95°F 13°C - 35°C	Storage Temperature Température de stockage Lagertemperatur Temperatura de almacenamiento Temperatura di immagazzinaggio Opslagtemperatuur Oppbevaringstemperatur Temperatura de conservação
REF	Catalog Number Numéro de référence Katalognummer Número de catálogo	Numero di catalogo Catalogusnummer Katalognummer Número de catálogo
	Manufacturer Fabricant Hersteller Fabricante	Produttore Fabrikant Produsent Fabricante
	Caution Attention Vorsicht Cuidado	Attenzione Let op Obs! Cuidado
	Consult Instructions for Use Consulter le mode d'emploi ! Gebrauchsanleitung konsultieren! Consulte las Instrucciones de uso Consultare le Istruzioni per l'uso Raadpleeg de gebruiksaanwijzing! Se bruksanvisningen! Consultar as Instruções de Utilização!	
LOT	Lot Number Numéro de lot Chargenbezeichnung Número de lote	Numero di lotto Partijnummer Partinummer Número de lote
	Use By Utiliser avant Verwendbar bis Fecha de caducidad	Data di scadenza Uiterste gebruiksdatum Brukes før Validade
STERILEEO	Sterilized using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Mit Ethylenoxid sterilisiert Esterilizado mediante óxido de etileno Sterilizzato con ossido di etilene Gesteriliseerd met ethyleenoxide Sterilisert med etylenoksid Esterilizado por óxido de etileno	
EC REP	Authorized Representative in the European Community Représentant agréé dans l'Union européenne Autorisierte EU-Vertretung Representante autorizado en la comunidad europea Rappresentante autorizzato per l'Unione Europea Erkend vertegenwoordiger in de Europese Gemeenschap Autorisert representant i EU Representante Autorizado na Comunidade Europeia	
	Do not use if the packaging is damaged. Ne pas utiliser si le conditionnement est endommagé. Bei beschädigter Verpackung nicht verwenden. No debe utilizarse si el envase está dañado Non utilizzare se la confezione è danneggiata. Niet gebruiken als verpakking beschadigd is. Må ikke brukes hvis emballasjen er skadet. Não utilizar se a embalagem estiver danificada.	
	Do Not Reuse Ne pas réutiliser Nicht wiederverwenden No reutilizar	Non riutilizzare Niet opnieuw gebruiken Skal ikke gjenbrukes Não reutilizar
	Do not resterilize Ne pas restériliser. Nicht resterilisieren. No lo vuelva a esterilizar.	Non risterilizzare. Niet hersteriliseren. Må ikke steriliseres igjen. Não reesterilizar.
Manufactured in Mexico for:		Manufactured in Mexico for: Fabriqué au Mexique pour : Hergestellt in Mexiko für: Hecho en México para: Prodotto in Messico per: Vervaardigd in Mexico voor: Produsert i Mexico for: Fabricado no México para:

Graft Preparation and Delivery Device

EN

Before using, read the following information:

Indications for Use:

The Graft Preparation and Delivery Device is intended for the delivery of allograft and autograft or synthetic bone graft materials to an orthopaedic surgical site. In addition, it is designed to facilitate the premixing of bone graft materials with fluids such as I.V. fluids, blood, plasma, platelet rich plasma, bone marrow or other specified blood components as deemed necessary by the clinical use requirements.

Precautions:

There are no specific warnings, precautions, or adverse effects associated with the use of this device.

- Instructions for Use:**
1. Thread barrel of 15±0.5ml graft chamber into graft chamber cap, located in tray. Place funnel in proximal end of graft chamber.
 2. Fill chamber with the Prescribed quantity of graft material (allograft, autograft, demineralized bone, or synthetic bone graft material). Remove Funnel.
 3. Insert plunger into graft chamber and depress until it meets or rests slightly above graft material.
 4. To hydrate graft material, remove red end cap, twist blending connector with dual syringe Applicator or a single syringe onto distal luer of graft chamber. Depress plunger(s) of syringes to transfer fluids to 15±0.5ml graft chamber.
 5. Place 15±0.5ml graft chamber vertically in tray until ready for use.
 - 6-7. To deliver graft material to surgical site, remove end cap and depress plunger until fully extended beyond graft chamber.
 8. Dispose of as biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

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