
BIOMET[®]

BIOLOGICS

RIGID ENDOSCOPIC AEROSOL APPLICATOR TIP

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning
Instruções de Utilização

DIST

Biomet Biologics
56 East Bell Drive
P.O. Box 587
Warsaw IN 46581, USA

Manufactured in Mexico for:

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977

CE 0297


EC REP

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




6

RIGID ENDOSCOPIC AEROSOL APPLICATOR TIP

Before using,
read the following information:

**** IMPORTANT NOTE ****

This device is intended for use with a syringe or compatible applicator with a standard luer lock and Regulator. These must be ordered separately. Make sure you have all parts before proceeding.

DESCRIPTION / INDICATIONS FOR USE:

The Rigid Endoscopic Aerosol Applicator Tip is intended for the application of two non-homogeneous liquids.

CONTRAINDICATIONS: None known.

PRECAUTIONS:

- Connect to regulated CO2 gas source; maximum input pressure not to exceed 100 psi (7 bar). Air or gas embolism has occurred related with the use of spray devices employing pressure regulated gas at higher than recommended pressures in close proximity to tissue surface. When using pressurized spray devices changes in blood pressure, pulse, oxygen, saturations, and end tidal CO2 should be monitored.
- If lever force is applied to the Applicator, grasp the upper portion of the stainless steel shaft. Do not use the syringe assembly as a hold point.
- Ensure spray tip is secured to prevent accidental tip detachment.
- Retain Tip Alignment Tool which holds a second Replacement Tip provided for use if the first tip becomes occluded.
- (Figure 6) If Replacement Tip becomes occluded during use, (indicated by no spray or no movement of the ball in the gas flow gauge) use the Tip Alignment Tool (which the installed Tip came from) to remove plugged Tip by unscrewing counterclockwise.

INSTRUCTIONS FOR USE:

1. [Note: Thoroughly read and refer to the Regulator IFU for proper set up and use.] (Figure 1/1.1) Turn on the Regulator gas supply and adjust gas source pressure to 100 psi (7 bar). The pressure gauge should indicate 18-22 psi (1.24-1.5 bar) or in the green zone. (Figure 2/2.1/2.2) Adjust the Regulator flow to 1.0-2.5 liters per minute. Check the gas flow by noting the height of the ball in the flow gauge while holding the foot switch down.
2. REMOVE the applicator, tubing, and tip from its sterile package using sterile technique.
3. (Figure 3) Use Tip Alignment Tool to attach the tip to the Applicator.
4. (Figure 4/4.1/4.2) Attach gas supply line (clear luer connector) to Applicator (4.1) and attach the patient vent line (red luer connector) to available female luer on trocar cannula vent valve (4.2).
 - The twin tube assembly can be easily separated as necessary to allow freedom of movement.
5. (Figure 5) Attach spray set to Regulator.
 - Attach the large clear fitting on tubing set to male luer connection. Attach the blue luer connection on tubing set to the female luer connection.
6. Assemble the filled syringes to the Applicator and install the syringe holder and plunger link per Applicator Assembly instructions.
7. Apply pressure to the foot switch to start gas flow prior to applying liquids and ensure flow by checking the gauge. If flow gauge does not move when the foot pedal is depressed, replace the tip.
8. While activating foot switch, dispense liquids by depressing the syringe plungers with VERY SLOW STEADY PRESSURE. Do not spray at a distance closer to the surface of tissues than 3cm at a maximum flow rate of 25 psi. To stop spray delivery, release pressure on syringe plungers while maintaining gas flow by holding foot down for 3-5 seconds to clear out the Applicator.

9. Dispose of as Biohazard waste. Device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness, or death of a patient.

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