BIOMET[®]

BIOLOGICS

Regulator

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning

DIST

BIOMET[®]

BIOLOGICS

56 EAST BELL DRIVE
P.O. BOX 587
WARSAW, IN 46581 USA

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul, MN 55121 USA
651 452-1977

CE 0297

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
Regulator Parts List	
1	Hose assembly
2	Foot switch
3	Pressure adjustment knob
4	Pressure gauge
5	Luer gas outlet

1. INTENDED USE

The Control Unit Regulator and Aerosol Applicators are intended for use in applying two fluids to a treatment site using compressed air or an inert gas.

2. CONTRAINDICATIONS

The Regulator is contraindicated for intravascular use and for procedures involving insufflating a confined body cavity.

3. WARNINGS/ PRECAUTIONS

The user should be familiar with the contents of these instructions prior to using the device. Only qualified personnel should operate this device. Caution must be used when applying product using pressurized gas.

- Air or gas embolism has occurred with the use of spray devices employing a pressure regulator to administer fibrin sealants. This event appears to be related to the use of the spray device at higher than recommended pressures and in close proximity to the tissue surface. When applying sealants using a spray device, be sure to use the flow rate recommended in the Instructions for use.
- To avoid possible gas embolism, do not spray directly into circulatory pathways. Any application of pressurized gas is associated with a potential risk of gas embolism, tissue rupture or gas entrapment with compression, which may be life threatening.

Be sure to take appropriate measures to address these risks by observing these recommendations:

- Do not spray at a distance closer to the surface of tissues than 3 cm at a maximum pressure of 25 psi.
- When using pressurized spray devices, changes in blood pressure, pulse, oxygen saturations, and end tidal CO₂ should be monitored because of the possibility of occurrence of gas embolism.
- Review Applicator instruction for use for additional warnings and precautions.
- Connect only to a pre-regulated supply source of medical gas (inert gas or air).
- Recommended supply pressure to Regulator should be set to 100 ± 5 psi ($6.89 \pm .34$ bar). Maximum supply pressure not to exceed

130 psi (8.96 bar).

- Do not expose to temperatures below -4°F (-20°C), above 122°F (50°C), or to relative humidity (RH) above 95%.
- Do not immerse or soak any part of the device.
- Review Section 7 for additional information regarding cleaning and other special instructions.

4. DEVICE DESCRIPTION

The Regulator is a reusable, non-sterile device that reduces incoming gas pressure to an adjustable output from 0 - 25 psi (0 - 1.72 bar). The pressure gauge has a reading of 0 - 30 psi (0 - 2.07 bar) with an accuracy of ± 2 psi (0.14 bar).

4.1 Contents of Carton

Regulator assembly, foot switch and hose assembly; pole clamp (optional).

NOTE: The Regulator comes standard equipped with a DISS (Female Hand Tight) medical air or CO₂ hose fitting, depending on Regulator model. This will mate to the DISS (Male) medical air or CO₂ fitting in the Operating Room.

NOTE: Operating Rooms may vary in supply fittings for providing pressurized gases. It may be necessary for the user to provide the appropriate adapter to connect the Regulator hose assembly to operating room supply fitting.

4.2 Set-up

Note: Allow Regulator to adjust to room temperature, if necessary.

After removing device from the package, inspect device for any mechanical or physical damage that could impair operation.

- Regulator: Attach Hose Assembly (1) to pre-regulated supply source and set supply pressure to 100 ± 5 psi (6.89 ± 0.34 bar).
- Foot switch: Place the foot switch (2) flat in a convenient location.
- Pole clamp installation (optional): Align the holes in the flat plate of the pole clamp with the threaded holes on the bottom face of the regulator. The clamp should be above the plate when in place.
- Thread the knobs (included with pole clamp) through each of the holes and tighten securely.

5. FUNCTIONAL CHECK

Conduct a functional check prior to first time use and minimally on a semi- annual basis.

5.1 Verify Regulator Pressure Output

- Unlock the pressure adjustment knob (3) by lifting it away from the surface of the Regulator.

- While depressing the foot switch, turn the pressure adjustment knob clockwise and counterclockwise and verify pressure output is adjustable from 0 to 25 psi (0 - 1.72 bar) as indicated on the pressure gauge (4).
- If pressure output is verified, the Regulator is ready for use (proceed to next section).

5.2 If the pressure output is not verified:

- Remove Regulator from service.
- Contact Micromedics Customer Service using the contact information on the device.

6. INSTRUCTIONS FOR USE

6.1 Attach the hose assembly to the gas source in the operating room (see Section 4.1).

6.2 Set the gas pressure on the Regulator:

- Unlock () the pressure adjustment knob by lifting it away from the surface of the regulator box.
- Rotate the pressure adjustment knob to the desired pressure setting within the green zone.
- Lock () the knob by pushing it toward the surface of the regulator box until it locks with a snap.

6.3 Attach the Applicator to the Regulator (5).

6.4 Control gas flow (on/off) to the Applicator by depressing or releasing the Foot Switch.

6.5 If the spray pattern for the fluids applied is not acceptable, increase the pressure setting up to 25 psi (1.72 bar) while monitoring the pressure gauge, Applicator spray pattern and Applicator fluid flow volume.

Coarse spray (large droplets) – low pressure (15 psi/1.03 bar)

Fine spray (small droplets) – high pressure (20 psi/1.38 bar)

Do not exceed 25 psi (1.72 bar)

7. CLEANING

The exterior surfaces of the Regulator, foot switch, and hose can be cleaned with many hospital approved (pH 7 neutral) disinfectants.

- In order for the Regulator to operate safely and in accordance with the intended use, the Regulator must be cleaned following each use and maintained free of dirt and other debris.
- CAUTION: Do not use high pH (acid) or low pH (caustic cleansers).
- CAUTION: Do not allow cleaning solution to enter outlet (5) fitting on front panel. Use Luer

caps or cover fittings with a lint free cloth.

- CAUTION: Do not allow cleaning solution to enter the Hose Assembly (1); protect with a lint free cloth.
- CAUTION: Do not allow disinfectants to pool anywhere on the device.
- CAUTION: Do not immerse or soak the Regulator in liquids.
- CAUTION: Do not use abrasive pads or cleansers.
- CAUTION: Do not process in automatic washer or steam sterilizer.

8. MAINTENANCE

NOTE: REGULATOR INTERNAL COMPONENTS ARE NOT USER SERVICEABLE.

- Perform a maintenance check annually.
- While depressing the foot switch, adjust the pressure indicated on gauge between 0 psi (0.0 bar) and 25 psi (1.72 bar).
If the gauge does not indicate down to 0 psi (0.0 bar) or moves erratically, remove the Regulator from service and contact Micromedics for further instructions.

9. STORAGE

Conditions:

- Ambient indoor 55° F (13° C) to 95° F (35° C), 0% to 95% RH continuous.
- Protect from damage, loose debris, dust, frost and direct sunlight.
- During first time use, after storage, and minimally on a semi-annual basis perform the Functional Check (See Section 5.0).

10. DISPOSAL

The Regulator has no electrical components or hazardous materials in its construction. Clean the Regulator to remove any potential biohazard contamination prior to disposal.

- Do not dispose of the Regulator as unsorted municipal waste.
- Collect the Regulator as used medical waste.
- Use collection and return systems available to you.

11. WARRANTY

The manufacturer warrants devices against all defects in material and workmanship that occur within 1 year after the date of sale. Warranty claims will be processed by the manufacturer under the following conditions:

- Inform Micromedics immediately of malfunction. Please provide the lot number

of the device when reporting.

- Observe all instructions regarding storage and return of the device.
- Present a legible copy of the invoice for the device in question which clearly shows date of sale.
- Describe in as much detail as possible suspected defects or malfunctions experienced by user.

ATTENTION: THE WARRANTY WILL NOT BE APPLICABLE IF WARNINGS, INTENDED USE, AND CLEANING INSTRUCTIONS FOR THE DEVICE ARE NOT FOLLOWED AS STATED HEREIN. THE MANUFACTURER AND THE SUPPLIER OF THE DEVICE DENY LIABILITY, IF:

- The device is not used in accordance with instructions for use.
- The defect was not caused by manufacturer or persons and service companies authorized by the same.

12. RETURN OF DEVICES

Reasons for return:

1. Results of the first-time use are unsatisfactory or unsafe.
2. Defective device:
 - Contact Micromedics Customer Service using the information provided on the IFU to obtain a Return Material Authorization (RMA) number prior to returning. Product returns are not accepted without an RMA number.
 - When returning a device, clean it carefully and pack it in the original package. If the original package is no longer available, request appropriate packaging through Micromedics Customer Service.

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