The Medicines Company

RaplixaReg™ Air Regulator

REF PF-FC-001

DIST

ProFibrix BV (The Medicines Company)
Darwinweg 24
2333 CR Leiden
The Netherlands

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul, MN 55121-1385, USA
651-452-1977
800-624-5662 USA
info@nordsonmicromedics.com
www.nordsonmicromedics.com

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



RaplixaReg™ Air Regulator

Before using, read the following information:

INTENDED USE

The RaplixaReg™ Air Regulator is intended for use with the RaplixaSpray device in applying biomaterials to a treatment site using compressed air or an inert gas.

CONTRAINDICATIONS

The RaplixaReg™ Air Regulator is contraindicated for intravascular use and for procedures involving insufflating a confined body cavity.

WARNING AND PRECAUTIONS

Connect only to compressed air or other inert gas. Do not use with oxygen, carbon dioxide, nitrogen or other flammable gases.

Supply pressure must not exceed 130 psi (9 bar).

Any application of pressurized gas may be associated with a potential risk of gas embolism, tissue rupture or gas entrapment with compression, which may be life-threatening. Events of air or gas embolism have occurred that were related to the use of other marketed spray devices administrating fibrin sealants at higher than recommended pressures and in close proximity to the tissue surface. Caution must be used to exclude these risks, by control of maximum pressure (1.7 bar - 25 psi) and minimum distance (5 cm - 2 inches). Be sure to use the pressure recommended in the Instruction for Use of the RaplixaSpray. Do not touch the tip against the bleeding site and do not spray into circulatory pathways. When using pressurized spray devices, changes in blood pressure, pulse, oxygen saturations, and end tidal CO₂ should be monitored.

ADVERSE REACTIONS

There are no known adverse reactions.

DISS MEDICAL AIR SUPPLY REQUIREMENTS

Note: The RaplixaReg™ Air Regulator comes standard equipped with an ISO 5359-compliant hose (black and white) and DISS (Female Hand Tight) Medical Air hose fitting. This will mate to DISS (Male) Medical Air fitting in the Operating Room.

Hose color requirements and supply fittings may vary by geography and/or operating rooms (i.e., the standard hose color in the USA is yellow) for providing pressurized Medical Air. It may be necessary for the user to change the hose and/or supply fitting to connect the Air Regulator Medical Air hose assembly to operating room fitting.

7 bar (100 psi) gas input pressure provides optimum performance for the RaplixaReg™ Air Regulator. Maximum gas supply pressure not to exceed 9 bar (130 psi).

OPERATING INSTRUCTIONS

Note: Allow the Air Regulator to adjust to room temperature if necessary.

Conditions: Ambient indoor +55°F (+13°C) to +95°F (+35°C), 0% to 60% RH.

The following steps are to be taken:

- 1 Attach the DISS Medical Air hose from the RaplixaReg™ Air Regulator to the Operating Room Medical Air supply (DISS Medical Air). (Not shown.)
- 2 Set the gas pressure on the Air Regulator.
2A Unlock (🔓) the pressure adjustment knob by lifting it away from the surface of the regulator box.
2B The knob should be rotated to the optimal pressure setting of 1.5 bar (22 psi).
2C Lock (🔒) the knob by pushing it toward the surface of the regulator box until it locks with a snap.
- 3 Attach the ProFibrix device to the regulator.
- 4 Start the flow of air to the ProFibrix device by turning the Air Regulator switch to the "ON" (I) position.
- 5 Turn the Air Regulator switch to "OFF" (O) position to stop flow of gas. (Not shown.)

A RaplixaReg™ Air Regulator pressure set in excess of 1.7 bar (25 psi) may cause safety relief valve to begin actuation. The safety relief valve will reset automatically when pressure is temporarily reduced to less than 0.5 bar (7 psi).

MAINTENANCE

On an annual basis, it is recommended to perform a maintenance check. With the Air Regulator turned on (I), adjust pressure indicated on gauge between 0 psi/0 bar and 22 psi/1.5 bar. If the gauge does not indicate down to 0 psi/0 bar or moves erratically, remove the Air Regulator from service and contact the supplier for further instructions.

CLEANING

The Air Regulator box and gas supply hose may be cleaned by wiping with soap and water. The materials used in the construction of the Air Regulator box and gas supply hose are unaffected by brief exposure to disinfectants used in the operating room.

STORAGE

Conditions for storage:

- Ambient indoor +55°F (+13°C) to +95°F (+35°C), 0% to 95% RH continuous.
- Protect from damage, loose debris, dust, frost and direct sunlight.

DISPOSAL

The RaplixaReg™ Air Regulator has no electrical components or hazardous materials in its construction. Clean the Air Regulator to remove any potential biohazard contamination prior to disposal.

- Do not dispose of the Air Regulator as unsorted municipal waste.
- Collect the Air Regulator as used medical waste.
- Use collection and return systems available to you.

WARRANTY

Products will be free of defects in material and workmanship for a period of two (2) years from the date of shipment. In the event a product is defective and notice of defect is given to Micromedics within two (2) years from the date of shipment, the Manufacturer will, at its option, repair or replace the defective product or credit the Buyer with the purchase price paid. The warranty is not applicable for damage or defect NOT caused by the manufacturer.

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