The Medicines Company

RaplixaReg™ CO2 Regulator

REF PF-FC-009

DIST

ProFibrix BV (The Medicines Company)
Darwinweg 24
2333 CR Leiden
The Netherlands

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul, MN 55121-1385, USA
651-452-1977
800-624-5662 USA
info@nordsonmicromedics.com
www.nordsonmicromedics.com

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



RaplixaReg™ CO2 Regulator

⚠ Before using, read the following information:

INTENDED USE

The RaplixaReg™ CO2 Regulator is intended for use with the RaplixaSpray device in applying biomaterials to a treatment site using compressed CO2 gas.

CONTRAINDICATIONS

The RaplixaReg™ CO2 Regulator is contraindicated for intravascular use and for procedures involving insufflating a confined body cavity.

WARNING AND PRECAUTIONS

Connect only to compressed CO2 gas. Do not use with compressed air, oxygen, nitrogen or other flammable gases.

Supply pressure must not exceed 130 psi (9 bar).

Any application of pressurized gas may be associated with a potential risk of gas embolism, tissue rupture or gas entrapment with compression, which may be life-threatening. Events of air gas embolism have occurred that were related to the use of other marketed spray devices administering fibrin sealants at higher than recommended pressures and in close proximity to the tissue surface. Caution must be used to exclude these risks, by control of maximum pressure (1.7 bar - 25 psi) and minimum distance (5 cm - 2 inches). Be sure to use the pressure recommended in the Instruction for Use of the RaplixaSpray. Do not touch the tip against the bleeding site and do not spray into circulatory pathways. When using pressurized spray devices, changes in blood pressure, pulse, oxygen saturations, and end tidal CO2 should be monitored.

ADVERSE REACTIONS

There are no known adverse reactions.

DISS CO2 SUPPLY REQUIREMENTS

Note: The RaplixaReg™ CO2 Regulator comes standard equipped with an ISO 5359-compliant hose (grey) and DISS (Female Hand Tight) CO2 hose fitting. This will mate to DISS (Male) CO2 fitting in the Operating Room.

Operating Rooms may vary in supply fittings for providing pressurized CO2. It may be necessary for the user to provide the appropriate adapter to connect CO2 Regulator CO2 hose assembly to operating room supply fitting.

7 bar (100 psi) gas input pressure provides optimum performance for the RaplixaReg™ CO2 Regulator. Maximum gas supply pressure not to exceed 9 bar (130 psi).

OPERATING INSTRUCTIONS

Note: Allow the CO2 Regulator to adjust to room temperature if necessary.

Conditions: Ambient indoor +55°F (+13°C) to +95°F (+35°C), 0% to 60% RH.

The following steps are to be taken:

- 1 Attach the DISS CO2 hose from the RaplixaReg™ CO2 Regulator to the Operating Room CO2 supply (DISS CO2). (Not shown.)
- 2 Set the gas pressure on the CO2 Regulator.
2A) Unlock (🔓) the pressure adjustment knob by lifting it away from the surface of the regulator box.
2B) The knob should be rotated to the optimal pressure setting of 1.5 bar (22 psi).
2C) Lock (🔒) the knob by pushing it toward the surface of the regulator box until it locks with a snap.
- 3 Attach the ProFibrix device to the regulator.
- 4 Start the flow of CO2 to the ProFibrix device by turning the CO2 Regulator switch to the "ON" (I) position.
- 5 Turn the CO2 Regulator switch to "OFF" (O) position to stop flow of gas. (Not shown.)

⚠ A RaplixaReg™ CO2 Regulator pressure set in excess of 1.7 bar (25 psi) may cause safety relief valve to begin actuation. The safety relief valve will reset automatically when pressure is temporarily reduced to less than 0.5 bar (7 psi).

MAINTENANCE

On an annual basis, it is recommended to perform a maintenance check. With the CO2 Regulator turned on (I), adjust pressure indicated on gauge between 0 psi/0 bar and 22 psi/1.5 bar. If the gauge does not indicate down to 0 psi/0 bar or moves erratically, remove the CO2 Regulator from service and contact the supplier for further instructions.

CLEANING

The CO2 Regulator box and gas supply hose may be cleaned by wiping with soap and water. The materials used in the construction of the CO2 Regulator box and gas supply hose are unaffected by brief exposure to disinfectants used in the operating room.

STORAGE

Conditions for storage:

- Ambient indoor +55° F (+13°C) to +95°F (+35°C), 0% to 95% RH continuous.
- Protect from damage, loose debris, dust, frost and direct sunlight.

DISPOSAL

The RaplixaReg™ CO2 Regulator has no electrical components or hazardous materials in its construction. Clean the CO2 Regulator to remove any potential biohazard contamination prior to disposal.

- Do not dispose of the CO2 Regulator as unsorted municipal waste.
- Collect the CO2 Regulator as used medical waste.
- Use collection and return systems available to you.

WARRANTY

Products will be free of defects in material and workmanship for a period of two (2) years from the date of shipment. In the event a product is defective and notice of defect is given to Micromedics within two (2) years from the date of shipment, the Manufacturer will, at its option, repair or replace the defective product or credit the Buyer with the purchase price paid. The warranty is not applicable for damage or defect NOT caused by the manufacturer.

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