
CO₂ Regulator

REF AX-610

Catalog Number

Caution

Consult Instructions
for Use

Do not use if package is damaged

DIST

Distributed By

Actamax Surgical Materials LLC
735 Pennsylvania Drive
Exton, PA 19341 USA
www.actamax.com

Manufacturer

Micromedics, Inc.
1270 Eagan Industrial Road
St. Paul MN 55121-1385, USA
651-452-1977
800-624-5662
www.nordsonmicromedics.com



CE 0297

EC REP

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1	DISS CO ₂ Gas Connection
2	Gas Outlet to Applicator
3	Patient Vent Connection
4	Flow Meter
5	Toggle Switch (Gas Flow ON/OFF)
6	Pole Clamp
7	Flow Control

CO₂ Regulator

REVIEW ACTAMAX™ ADHESION BARRIER INSTRUCTIONS FOR USE FOR ADDITIONAL INFORMATION, WARNINGS AND PRECAUTIONS.

1.0 INTENDED USE

The AX-610 CO₂ Regulator is a re-usable, nonsterile device intended for regulating gas flow of medical grade vaporized carbon dioxide to assist in the application of biomaterials.

2.0 CONTRAINDICATIONS

None.

3.0 WARNINGS/PRECAUTIONS

The user should be familiar with the contents of these instructions prior to using the device.

Only qualified personnel should operate this device.

Do not use this device with any spray applicator other than spray applicators recommended by Actamax for use with the Actamax™ Adhesion Barrier.

Caution must be used when applying product using pressurized gas.

- Air or gas embolism has occurred with the use of spray devices employing a pressure regulator to administer fibrin sealants. This event appears to be related to the use of the spray device at higher than recommended pressures and in close proximity to the tissue surface. When applying biomaterials using a spray device, be sure to verify the flow rate recommended in the Instruction for Use.
- To avoid possible gas embolism, do not spray directly into circulatory pathways. Any application of pressurized gas is associated with a potential risk of air embolism, tissue rupture or gas entrapment with compression, which may be life threatening.

Be sure to take appropriate measures to address these risks by observing these recommendations:

- Do not spray at a distance closer to the surface of tissues than 3 cm at a maximum flow rate of 4.5 liters per minute.
- When using pressurized spray devices, changes in blood pressure, pulse, oxygen saturations, and end tidal CO₂ should be monitored because of the possibility of occurrence of gas embolism.

Connect to regulated CO₂ gas source only.

Supply pressure to device shall be set to 100 ± 5 psi (6.89 ± 0.35 bar) and 10 liters per minute minimum flow.

Do not connect to unregulated supply pressure greater than 250 psi (17.2 bar).

Do not expose to temperatures below -4°F (-20°C), above 122°F (50°C), or to relative humidity (RH) above 95%.

Do not immerse or soak any part of the device.

Review Section 8.0 for additional information regarding cleaning, maintenance and other important information for use.

4.0 DESCRIPTION OF OPERATION

The AX-610 CO₂ Regulator is a reusable, nonsterile device consisting primarily of a gas pressure regulator, flow control valve, flow meter and an ON/OFF switch. The internal pressure regulator is a fixed-output device that reduces 100 psi ± 5 (6.89 bar ± 0.35) incoming supply gas pressure to the device to a preset 19 psi ± 2 (1.38 bar ± 0.13) output pressure. The flow output is adjustable. A safety relief valve inside the Regulator will safely vent inside the housing at pressures above 22 psi (1.5 bar) with an audible signal. The pressure relief valve will reset automatically when the over pressure condition is removed.

4.1 CONTENTS OF CARTON

Regulator assembly (including pre-mounted pole clamp), Gas Hose Assembly with DISS Female CO₂ Connection.

5.0 SET UP

Conditions: 55° F (13° C) to 95° F (35° C), 0% to 60% RH

Note: Allow Regulator to adjust to room temperature before use.

The AX-610 CO₂ Regulator device comes with a pre-mounted Pole Clamp (6) that can be used to mount the device on a stable pole.

After removing device from the package, inspect device for any mechanical or physical damage that could impair operation. Inspect all labels and/or markings for legibility and completeness.

5.1 SET UP REGULATOR

- Attach Gas Hose Assembly (DISS Female CO₂ Connection) (1) to CO₂ gas source and set gas supply pressure to 100 psi ± 5 (6.89 bar ± 0.35) and 10 liters per minute minimum flow.
- To turn on gas flow, rotate the ON/OFF switch (5) clockwise.
- To turn off gas flow, rotate the ON/OFF switch (5) counter-clockwise.

6.0 FUNCTIONAL CHECK

During first time use, after storage and minimally on a semiannual basis, conduct a functional check.

6.1 VERIFY REGULATOR FLOW

- Turn on the ON/OFF switch (5) and verify that the Flow Meter (4) reads flow through the regulator.
- If the Flow Meter (4) does not read flow through the regulator, verify that the Flow Control Valve (7) is open by rotating it in a counter-clockwise direction.

- Turn off the ON/OFF switch (5) and verify flow control ball does not stick at top of Flow Meter (4).

Note: After performing above procedure, the AX-610 CO₂ Regulator is ready for use.

If the flow is not verified:

- Immediately remove the AX-610 CO₂ Regulator from service.
- Contact Actamax using the contact information on the rear panel of the device.

6.2 SAFETY RELIEF VALVE OPERATION

CAUTION: Pressure above 22 psi / 1.5 bar will cause the Regulator's internal safety relief valve to release excess gas pressure inside Regulator assembly with an audible signal.

Should the safety relief valve actuate at anytime:

- Confirm incoming gas supply pressure to the device is set correctly to 100 psi $\pm$ 5 (6.89 bar $\pm$ 0.35)
- Should the safety relief valve continue to actuate after the device input pressure is confirmed to be set correctly:
- Immediately remove AX-610 CO₂ Regulator from service.
- Contact Actamax using the contact information on the rear panel of the device.

7.0 STORAGE

Conditions:

- Ambient indoor +55° F (+13° C) to +95° F (+35° C), 0% to 95% RH continuous.
- Protect from damage, loose debris, dust, frost and direct sunlight.

8.0 CLEANING AND SPECIAL INSTRUCTIONS

8.1 CLEANING

- The AX-610 CO₂ Regulator housing, and hose can be cleaned with many hospital approved (pH 7 neutral) disinfectants
- In order for the device to operate safely and in accordance with the intended use, the Regulator must be cleaned following each use and maintained free of dirt and other debris.

8.2 SPECIAL INSTRUCTIONS FOR CLEANING

- CAUTION: Do not use high pH (caustic) or low pH (acid) cleansers.
- CAUTION: Do not allow cleaning solution to enter Patient Vent (3) or Gas Outlet (2) fittings on front panel. Use Luer caps or cover fittings with a lint free cloth.
- CAUTION: Do not allow cleaning solution to enter the Gas Hose Assembly (DISS Female CO₂ connection) (1); cover with a lint free cloth.

- CAUTION: Do not allow disinfectants to pool anywhere on the device such as on Flow Meter (4).
- CAUTION: Do not immerse or soak the entire unit.
- CAUTION: Do not use abrasive pads or cleansers.
- CAUTION: Do not process in automatic washer or steam sterilizer.

9.0 SCHEDULED MAINTENANCE CHECK

ACTAMAX AX-610 CO₂ REGULATOR INTERNAL COMPONENTS ARE NOT USER SERVICEABLE.

9.1 MAINTENANCE SET UP

- Attach Gas Hose Assembly (DISS Female CO₂ Connection) (1) to CO₂ gas source and set gas supply pressure to 100 psi $\pm$ 5 (6.89 bar $\pm$ 0.35).
- Connect a calibrated 0 – 30 psi (0 - 2.1 bar) pressure gauge to the Gas Outlet (2).
- Turn on the ON/OFF switch (5) and verify calibrated pressure gauge reading agrees to 19 psi $\pm$ 2 (1.31 bar $\pm$ 0.13). Flow gauge should indicate 0 liters per minute after a momentary jump.
- Remove calibrated pressure gauge.

9.2 VERIFY PATIENT VENT

- Connect the tubing set from the spray applicator recommended by Actamax for use with Actamax™ Adhesion Barrier to Gas Outlet (2) and Patient Vent (3).
- Note: A short single lumen tube set with appropriate Male/Female Luer connections may be substituted for the above.
- Turn ON/OFF switch (5) on and verify that the flow meter indicates flow.
- Turn off ON/OFF switch (5), gas flow must stop and Flow Meter (4) goes to zero.
- Remove Tubing set from Gas Outlet (2) and Patient Vent (3).

9.3 CHECK FOR LEAKS

- Block Gas Outlet (2) with a luer cap, and check for internal leaks by observing the movement of Flow Meter (4) ball. Turn on the ON/OFF switch (5). The ball will remain at zero (after a momentary jump) if free of leaks.
- Remove Luer Cap from Gas Outlet (2).

10.0 DISPOSAL OF DEVICES

The AX-610 CO₂ Regulator has no electrical components or hazardous materials in its construction.

Clean the AX-610 CO₂ Regulator to remove any potential biohazard contamination prior to disposal.

Do not dispose of AX-610 CO₂ Regulator as unsorted municipal waste.

Collect AX-610 CO₂ Regulator separately as used medical waste.



Use collection and return systems available to you.

For more information on recovery or recycling of AX-610 CO₂ Regulator, contact Actamax using the contact information on the rear panel of the device.

11.0 TECHNICAL DATA

Outer dimensions:	L 216mm(8.5") x W216mm(8.5") x H165mm(6.5")
Weight:	2 kg (4.4 lbs)
Input Pressure:	100 psi ± 5 (6.89 bar ± 0.35)
Operating Conditions:	55° F (13° C) to +95° F (35° C) 0% to 60% RH
Flowmeter Accuracy:	If flowmeter reads between 3.5 - 4.5 LPM, actual flow is also within 3.5 - 4.5 LPM when connected to spray applicators recommended by Actamax for use with the Actamax™ Adhesion Barrier.

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