



QL® INFLATION Devices

Atrion, a Nordson MEDICAL Company, offers a versatile portfolio of high-pressure QL® Inflation Devices and Fluid Dispensing Syringes for a broad range of clinical applications including balloon catheter dilation, stent deployment, and high pressure fluid delivery.

These devices are available in multiple configurations to support diverse procedural needs. All products are manufactured in an ISO 13485, FDA registered facility and can be customized or private labeled to meet specific user requirements.

FEATURES & BENEFITS

- **Smart Design:**

- Ergonomic handle provides stability and precise control during use
- Quick-latch mechanism allows rapid deflation, even from maximum pressure
- Transparent barrel ensures clear visibility of contents and volume
- Design supports thorough and confident purging by medical personnel

- **Custom-Engineered Solutions:**

- Collaborative engineering support for specialized inflation devices, locking syringes, and fluid delivery systems
- Custom functionality, materials, and performance options available
- Flexible design and manufacturing to meet precise customer specifications

- **Custom Branding Options:**

- Available in multiple configurations to meet a wide range of clinical requirements

SERVICES & APPLICATIONS

- **Customer-Focused Services:**

- Built to your exact specifications
- Delivered on time
- Precision and reliability from production to delivery

- **Broad Applications**


- Cardiology
- Peripheral Vascular
- Structural Heart
- Urology
- Gastroenterology
- Neurology
- Orthopedics
- Cell Delivery
- Pulmonary
- ENT

SPECIALTY DEVICES FOR BALLOON CATHETER DILATIONS

QL1455

Ultra High-Pressure


14cc/55ATM



QL2540

Ultra High-Pressure

25cc/40ATM



QL4515

Reduced Deployment Force


45cc/15ATM



QL6015

Large Capacity

60cc/15ATM



QL6020PB

Large Capacity

60cc/20ATM



STOCK CONFIGURATIONS

Part Number	Size	Pressure	Type
96500	10cc	15 ATM	Quick Latch
96501	10cc	30 ATM	Quick Latch
96351	14cc	55 ATM	Quick Latch
96502	20cc	30 ATM	Quick Latch
96450	20cc	30 ATM	Press-to-Release
96358	25cc	30 ATM	Quick Latch
96311	25cc	40 ATM	Quick Latch
96313	40cc	15 ATM	Quick Latch
96518	45cc	15 ATM	Quick Latch
96314	60cc	15 ATM	Quick Latch
96522	60cc	20 ATM	Push Button